

Fast multidimensional NMR spectroscopy of proteins in solution

Development of new methods and application to
the study of structure and kinetics

Paul Schanda

Directeur de thèse : Bernhard Brutscher




*Laboratoire de RMN
Institut de Biologie Structurale
Grenoble*



Major challenges in structural biology

Structures




X-ray diffraction
NMR (20% of PDB)

Static

β -catenin in complex with the intrinsically
unfolded protein tcf4

Major challenges in structural biology

Structures




X-ray diffraction
NMR (20% of PDB)

Static

β-catenin in complex with the intrinsically unfolded protein tcf4

Dynamics



Spectroscopies

IR

UV/Vis

NMR

...

Scattering

Diffraction

Transient Interactions

Folding

Unfolded proteins

Coupled binding/folding

Catalysis

Misfolding and fibrils

...


**Ensembles of structures
changing with time**

Introduction

Nuclear spins report on the environment at hundreds of sites within the molecule

Each individual NMR-active nucleus (^1H , ^2H , ^{13}C , ^{15}N ,...) gives local information


- Chemical environment
- Structure
- Dynamics



Introduction

Requirements to observe individual atomic sites

1. Sensitivity



Only 1 out of ~10000 molecules contributes to NMR signal


Solutions

High B_0 -fields




400 → 900MHz
factor ≈ 3.0

Cryogenic probes



factor ≈ 2-4

Optimized pulse sequences





Sensitivity has advanced a lot in recent years making possible a significant reduction of experimental durations

Introduction

Requirements to observe individual atomic sites

2. Resolution



76 amino acids


Observation of signals of individual nuclei is hindered by signal overlap

Solution


Multidimensional NMR

Introduction

Multidimensional NMR



... correlates the frequencies of several interacting nuclear spins and spreads the signals over a multidimensional space





Introduction

Time requirements of nD NMR



1D: seconds





2D: minutes




3D: hours



4D: weeks



Each point corresponds to one signal acquisition ($\approx 1-2$ seconds)



In many cases, experimental durations are dictated by the sampling requirement, not by the sensitivity.

Limitations of nD NMR related to long experimental durations

- * **Sample stability** and long acquisition times ?
- * **High throughput** capabilities for structural genomics projects ?
- * **High-dimensional data sets** as required for proteins with strong overlap (unfolded proteins) require unrealistically long times
- * **Spectral changes during a reaction** occurring on a seconds time scale cannot be followed with atomic resolution

Outline

Part 1: Methods for fast multidimensional NMR

Concepts

Sensitive fast-pulsing experiments: SOFAST HMQC, BEST experiments

Part 2: Applications


Application to the study of protein folding and unfolding

HET-SOFAST NMR: a fast tool for characterizing protein structure

Accelerating nD NMR

Concepts


1. Reduce the number of mandatory scans



Accelerating nD NMR

Concepts

1. Reduce the number of mandatory scans



Alternative ways of data sampling:

Non-linear data sampling of time domain

Hadamard NMR spectroscopy (frequency domain)


"Single scan" NMR (spatial/orientational encoding)

...

Accelerating nD NMR

Concepts

1. Reduce the number of mandatory scans



Alternative ways of data sampling:


Non-linear data sampling of time domain

Hadamard NMR spectroscopy (frequency domain)

"Single scan" NMR (spatial/orientational encoding)

...


2. Reduce the duration of each scan



Accelerating nD NMR

Concepts

1. Reduce the number of mandatory scans



Alternative ways of data sampling:


Non-linear data sampling of time domain

Hadamard NMR spectroscopy (frequency domain)





"Single scan" NMR (spatial/orientational encoding)

...





2. Reduce the duration of each scan







Considerations about sensitivity and interscan delay






Considerations about sensitivity and interscan delay



Considerations about sensitivity and interscan delay




Considerations about sensitivity and interscan delay



Accelerated longitudinal relaxation allows fast repetition of scans while retaining high sensitivity

Proton longitudinal relaxation in proteins





Most NMR experiments in liquids excite and detect **proton spins**.
The ^1H longitudinal relaxation determines the sensitivity.

Longitudinal relaxation in a two spin system


$$-\frac{d}{dt} \begin{pmatrix} I_{1z} - I_{1z}^0 \\ I_{2z} - I_{2z}^0 \end{pmatrix} = \begin{pmatrix} \rho & \sigma \\ \sigma & \rho \end{pmatrix} \begin{pmatrix} I_{1z} - I_{1z}^0 \\ I_{2z} - I_{2z}^0 \end{pmatrix}$$

ρ Auto-relaxation rate constant

σ Cross-relaxation rate constant



Proton longitudinal relaxation in proteins







The relaxation in a multi-spin system can be considered as a sum of two-spin interactions.

$$-\frac{d}{dt} \begin{pmatrix} I_{1z} - I_{1z}^0 \\ I_{2z} - I_{2z}^0 \\ \dots \\ I_{nz} - I_{nz}^0 \end{pmatrix} = \begin{pmatrix} \sum_j \rho_{1j} & \sigma_{12} & \sigma_{13} & \dots & \sigma_{1n} \\ \sigma_{21} & \sum_j \rho_{2j} & \sigma_{23} & \dots & \sigma_{2n} \\ \dots & \dots & \dots & \dots & \dots \\ \sigma_{n1} & \sigma_{n2} & \sigma_{n3} & \dots & \sum_j \rho_{nj} \end{pmatrix} \begin{pmatrix} I_{1z} - I_{1z}^0 \\ I_{2z} - I_{2z}^0 \\ \dots \\ I_{nz} - I_{nz}^0 \end{pmatrix}$$

- ρ Auto-relaxation rate constant
- σ Cross-relaxation rate constant

Amide ^1H relaxation after non-selective/selective excitation



Outline

Part 1: Methods for fast multidimensional NMR

Concepts

Sensitive fast-pulsing experiments: SOFAST HMQC, BEST experiments



Part 2: Applications

Application to the study of protein folding and unfolding

HET-SOFAST NMR: a fast tool for characterizing protein structure

2D ^1H - ^{15}N correlation experiment: **SOFAST-HMQC**

Band-**S**elective **O**ptimized-
Flip-Angle **S**hort-**T**ransient




Properties:

- * Selective pulses
- * Small number of rf pulses: reduces loss due to B_1 inhomogeneity
- * Good water suppression in only one scan


^1H - ^{15}N correlation experiment: **SOFAST-HMQC**

Band-**S**elective **O**ptimized-
Flip-**A**ngle **S**hort-**T**ransient




Properties:

- * Selective pulses
- * Small number of rf pulses: reduces loss due to B_1 inhomogeneity
- * Good water suppression in only one scan




^1H - ^{15}N correlation experiment: **SOFAST-HMQC**

Band-**S**elective **O**ptimized-
Flip-**A**ngle **S**hort-**T**ransient





Properties:

- * Selective pulses
- * Small number of rf pulses: reduces loss due to B_1 inhomogeneity
- * Good water suppression in only one scan
- * **Partial excitation (Ernst angle)**



Fast-pulsing regime




¹⁵N-ubiquitin (8.6 kDa, 25°C)

800 MHz, coldprobe


concentration: ~ 0.2 mM

Exp. time: 4 seconds



140°- SOFAST


Optimal-sensitivity regime



$t_{rec} \approx 0.2-0.3$ sec

SOFAST-HMQC performs best in terms of absolute sensitivity

α -lactalbumin (14 kDa, 25°C)
800 MHz, coldprobe
no ^{15}N enrichment, ~ 4.0 mM
(~ 10 μM ^{15}N -labeled protein)
Exp. time: 1 hour




120°- SOFAST

Can we go even faster ?

Combination with alternative encoding schemes


standard ^{15}N labeling



Can we go even faster ?


Combination with alternative encoding schemes

standard ^{15}N labeling




$N \approx 100$


Hadamard ^{15}N labeling




Spatial ^{15}N labeling



Hadamard encoded SOFAST HMQC



Spatially-encoded SOFAST-HMQC






Total acquisition time:
(8 interleaved scans)


1.4 seconds

3.4 mM ubiquitin
25°C, pH 6.0
800 MHz,
standard probe


Extension to 3D triple resonance experiments: **B**and-selective **E**xcitation **S**hort **T**ransient Experiments Example: 3D HNCA



3D HNCA experiment for resonance assignment in a few minutes



H-CA strips extracted from a 3D-BEST HNCA spectrum recorded in **15 min.**



Outline

Part 1: Methods for fast multidimensional NMR

Concepts



Sensitive fast-pulsing experiments: SOFAST HMQC, BEST experiments

Part 2: Applications



Application to the study of protein folding and unfolding

HET-SOFAST NMR: a fast tool for characterizing protein structure

Dynamics in proteins studied with atomic resolution by NMR



Dynamics in proteins studied with atomic resolution by NMR




Protein (un)folding, H/D exchange,...

SOFAST real-time 2D NMR

Required tools

Fast data acquisition

SOFAST HMQC




SOFAST real-time 2D NMR


Required tools

Fast data acquisition


SOFAST HMQC



Initiation of the kinetic event




Fast mixing inside the spectrometer within ca. 100 ms



Application I


Folding of α -lactalbumin followed in real time




α -lactalbumin 14 kDa
Native state pH 8

Application I

Folding of apo α -lactalbumin followed in real time



"Molten Globule"
pH 2




α -lactalbumin 14 kDa
Native state pH 8

Application I

Folding of apo α -lactalbumin followed in real time




2D snapshots of the folding process every 10 seconds



Application I




Folding of apo α -lactalbumin followed in real time

2D snapshots of the folding process every 10 seconds



Application I

Folding of apo α -lactalbumin followed in real time



No intermediate peaks observed


Same rate of disappearance of molten globule state and appearance of native state

Conclusion:

A single transition state ensemble controls folding


more folding...

Folding of an amyloidogenic protein





Work in progress

SOFAST-HMQC spectra recorded during refolding every 15 sec



more folding...

Folding of an amyloidogenic protein




Multiexponential intensity changes

Site-specific analysis of kinetic behavior will reveal nature of these states.


Application II

(Un)folding under native conditions studied by H/D exchange



Application II


(Un)folding under native conditions studied by H/D exchange




H/D exchange gives access to the "invisible" manifold of partially unfolded conformations under native conditions

Application II

Fast H/D exchange gives insight into unfolding





At **high pH** the observed exchange rates directly reflect (partial) unfolding rates



Application II

Unfolding of ubiquitin studied by EX1 H/D exchange




Global unfolding rate determined from fluorescence Chevron plots:

$$k_{NU} = 0.0084 \text{ s}^{-1} \quad \tau_{NU} = 119 \text{ s}$$

Bofill, R., ..., Searle, M.,
J. Mol. Biol. (2005) 349, 205


SOFAST HMQC-based H/D exchange experiments on 0.7mM ubiquitin at pH 12, 298 K



Application II

Unfolding of ubiquitin studied by EX1 H/D exchange

EX1 H/D exchange rates
plotted on ubiquitin's structure




- * Unfolding in ubiquitin is a heterogeneous process
- * Partial unfolding of the C-terminal face of the β -strand is faster
- * Gradual decrease of unfolding rates
- * The globular unfolding observed by fluorescence appears to be the unfolding of the N-terminal "folding nucleus"


Application II

Heterogeneous unfolding behavior of ubiquitin

EX1 H/D exchange



Cold denaturation



Ubiquitin in reverse micelles, -20°C
Babu et al. (2004) Nat Struct Mol Biol 11:352

Biological implication ?

Most thermostable H-bond: I3 → L15

Cordier and Grzesiek (2002), J.Mol.Biol, 715: 739.

Engineering a disulfide bridge between
4 and 66: 70-80% decrease of activity
in signaling proteolysis

Ecker et al. (1989) J Biol Chem 264, 1887

Outline

Part 1: Methods for fast multidimensional NMR

Concepts

Sensitive fast-pulsing experiments: SOFAST HMQC, BEST experiments

Part 2: Applications


Application to the study of protein folding and unfolding

HET-SOFAST NMR: a fast tool for characterizing protein structure

HET-SOFAST NMR

Amide proton relaxation exploited for structural characterisation

Towards a quantitative measure of structure:
what distinguishes **structured** and **unstructured**
polypeptide chains ?




- high proton density
- low water accessibility


- low proton density
- high water accessibility

HET-SOFAST NMR

Amide proton relaxation exploited for structural characterisation



non-selective(●)
selective(■)



$$\lambda = |\text{sat}/\text{ref}|$$




Amide proton relaxation depends on the structural context:

- * number of aliphatic protons close in space
- * local dynamics
- * water exchange of H^{N}

HET-SOFAST NMR

HET(erogeneity) SOFAST NMR *

Quantification of **proton density** and **water accessibility**
along polypeptide chains




Ubiquitin (2 mM, 25°C, 600 MHz)

Acquisition time: 10 s / spectrum

* Schanda, Forge & Brutscher,
Magn. Reson. in Chem. (2006) 44, S177.

HET-SOFAST NMR


2D HET-SOFAST NMR: characterizing structural compactness for assigned protein samples



Similar information to heteronuclear NOE
but obtained in significantly reduced experimental time



HET-SOFAST NMR

Structural and dynamic heterogeneities along polypeptide chains
prior to resonance assignment



HET-SOFAST NMR


Application: N-tailC, an unstructured protein with a small helical propensity



HET-SOFAST NMR allows to detect even partial structural preferences

HET-SOFAST NMR



Following subtle structural changes upon metal binding




α -lactalbumin

The structure of apo and holo forms are very similar
Metal binding was reported to induce some "structure rigidification"

HET-SOFAST NMR experiments can be used to follow small changes of protein "compactness"



CONCLUSIONS



Longitudinal relaxation optimized experiments (SOFAST, BEST) allow significant acceleration of nD NMR data acquisition

These approaches are fully compatible with alternative sampling methods (Hadamard, spatial encoding, reduced dimensionality,...)

Kinetic reactions in proteins can be probed simultaneously residue by residue on a time scale of seconds

^1H longitudinal relaxation can be used as a very sensitive probe for local structure and dynamics (HET-SOFAST NMR)

Merci

Bernhard Brutscher

Ewen Lescop

Mirjam Falge

Beate Bersch

Isabel Ayala

Cécile Giustini

Jean-Pierre Simorre

Adrien Favier

Catherine, Dominique, Fito,

Florence, Guillaume, Hélène,

Jérôme, Julien, Klaartje,

Laurence, Lionel, Malene, Martin,

Matthew, Monika, Nico, Pau,

Pierre, Phin, Rémy, Thomas





Thank you

Folding

Vincent Forge

Pulses, technical things

Eriks Kupce

Peter Sandor

« Ultrafast NMR »


Lucio Frydman

Maayan Gal


Longitudinal ^1H relaxation enhancement

Magnetic field dependence

simulation




experiment




Longitudinal relaxation optimized (selective) experiments become important especially for high field strength

Amide ^1H relaxation in proteins: Insight from simulation

Size of the molecule (tumbling rate)






Magnetic field strength



Fast pulsing


Aliphatic ^1H polarization recovery



Choice of selective pulses



Broadband Inversion Pulses

(Shaka & co, JMR 2001, 151, 269)






Time-optimized INEPT transfer

using PC9 (Kupce & Freeman, JMR 1993, 102A, 122)




Planar mixing without general rotation pulses


using EBURP2 (Geen & Freeman, JMR 1991, 93, 93)



Hadamard encoded SOFAST HMQC



Hadamard encoded SOFAST HMQC



Retrieving the 2D information
is achieved in only 2 scans !